

Mutector™ II

BRAF

V600 Mutation Detection Kit

Cat No. GP02-EK

User Manual V1.6

For research use only, not for use in diagnostic procedures

Limited Product Warranty

It is imperative that the users strictly adhere to this manual. Failure to do so will void TrimGen's guarantee of this product. TrimGen Corporation makes no other warranties of any kind, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose.

Notice to Purchaser

The Mutector™ II kit is provided as research use only, not for use in diagnostic procedures. The purchaser must determine the suitability of the product for their particular use.

The purchase of the Mutector™ II kit includes a limited, nonexclusive license to use the kit. This license does not grant rights to reproduce or modify the Mutector™ II kit for resale, or to use the Mutector™ II kit to manufacture commercial products without written approval of TrimGen Corporation. No other license, expressed, implied or by estoppels is granted.

Product Safety and Liabilities

When working with the kit reagents, always wear a lab coat, disposable gloves and protective goggles. TrimGen Corporation shall not be liable for any direct, indirect, consequential or incidental damages arising out of the misuse, the results of use, or the inability to use this product.

CONTENTS

Introduction	4
Overview of Mutector™ II Assay	5
Materials Provided	6
Materials Required	7
Equipment Required	7
DNA Sample Preparation	8
Sequencer Setup	8
Thermal Cycling Programs	9
Mutector™ II Assay Protocol	10
Data Analysis	15
Troubleshooting	16

Storage

Upon receipt of the kit, store at -20°C until use. At this temperature the reagents are stable for 8 months.

After first use, store all of the reagents at $2-8^{\circ}\text{C}$ and keep them protected from direct light. At this condition the reagents are stable for 2 months.

Introduction

The Mutector™ II BRAF V600 mutation detection kit (GP02-EK) uses TrimGen's proprietary technology called Shifted Termination Assay (STA) to detect mutations in codon V600 of the BRAF gene. This single tube assay detects the following mutations:

V600E (GTG>GAG)

V600A (GTG>GCG)

V600G (GTG>GGG)

V600K (GTG>AAG)


V600R (GTG>AGG)

The STA technology enriches mutation signals and is able to detect low levels of somatic mutations.

Mutector™ II has sequencing-like accuracy through multiple steps: (1) PCR amplification, (2) Selective sequence hybridization, (3) Sequence dependent termination and specific primer extension by STA.

* Shifted Termination Assay (STA)

Shifted Termination Assay (STA) technology is a proprietary multi-base primer extension method. The specially formulated ST reagents extend detection primers by multiple bases to increase signal strength and fragment size, creating mutation peaks that are easily distinguished from wild type. The STA technology has sequencing like accuracy and the mutation is confirmed by both peak color and fragment size. The method enriches mutation signal and detects low-level mutations missed by sequencing.



Overview of Mutector™ II Assay

Single tube assay for each sample

Step 1



PCR amplification

1-2 hours*

** Time varies depending on the type of thermal cyclers used*



Step 2



C-UP treatment

(PCR product clean-up)

30 min



Step 3



ST reaction

Mutation enrichment and detection

45 min*


** Time varies depending on the type of thermal cyclers used*



Load to sequencer



Step 4



Capillary Electrophoresis

Fragment analysis

30 min

Materials Provided:

The Mutector™ II *BRAF* Mutation kit (GP02-EK) contains reagents for 32 reactions.

Materials	Quantity
Master Mix	650 µL
PCR-P	50 µL
C-UP1	40 µL
C-UP2	40 µL
C-UP Buffer	350 µL
ST-BF*	400 µL
DP-BF	80 µL
CTL-BF	50 µL
Loading Buffer*	650 µL
TF - Filters	32
Collection Tubes	32 x 2

* **Light sensitive:** Keep these reagents protected from direct light.

Reagents Description:

Master Mix

Reagents for DNA amplification.

PCR-P

PCR primer mix for amplification of *BRAF* gene.

C-UP1 and C-UP2

Enzyme mix for cleanup of PCR products.

C-UP Buffer

Buffer for the C-UP enzymes.

ST-BF (Light sensitive)

Pre-mixed reagents for mutation enrichment and detection.

DP-BF

Pre-mixed detection primers.

CTL-BF

Pre-mixed mutant and wild type DNA. The controls are sufficient for 16 test runs.

Loading Buffer (Light sensitive)

Loading buffer for ABI capillary type sequencers and special fluorescence-labeled size standards.

Materials required:

0.2 ml PCR tubes (8-well strip tube)

DS-32 Matrix Standard Kit for the 3100 and 3130 Series Systems (one time set up. Applied Biosystems Cat. No. 4345831)

Equipment required:**Thermal Cycler:**

Any type of thermal cycler with a 0.2 ml tube block is acceptable for performing the Mutector™ II assay.

Sequencer:

Applied Biosystems capillary type Genetic/DNA Analyzer

Analysis Software:

Data Collection® software for ABI capillary sequencer
GeneMapper® for fragment analysis or GeneScan®

DNA Sample Preparation:

Any commercially available DNA extraction kit is acceptable.

Paraffin (FFPE) tissue samples

A high efficiency DNA extraction kit for FFPE sample is available from TrimGen.

Product information:

WaxFree™ DNA for 50 samples (Cat. # WF-50)

WaxFree™ DNA for 100 samples (Cat. # WF-100)

DNA concentration:

When using a column or bead DNA extraction method, adjust the final concentration of extracted DNA to **20-80 ng/ μ L**.

When using TrimGen's DNA preparation kit, follow the kit protocol to perform the PCR reaction.

Sequencer setup:

First time users should set up the analysis program for the sequencer (one time setup). After setup, the program can apply to all Mutector™ II tests for data analysis. Please choose either GeneMapper® or GeneScan® to analyze your data.

GeneMapper® Analysis

Step I. GeneMapper® Setup

www.trimgen.com/docs/PartI-GeneMapper-Setup.pdf

Step II. Data Collection® Software Setup

www.trimgen.com/docs/PartII-Data-Collection-Setup.pdf

Step III. Data Analysis Using GeneMapper®

www.trimgen.com/docs/PartIII-Data-Analysis-GeneMapper.pdf

GeneScan® Analysis

Step I. Data Collection® Software Setup

www.trimgen.com/docs/PartII-Data-Collection-Setup.pdf

Step II. GeneScan® Setup and Data Analysis

www.trimgen.com/docs/PartIV-Genescan.pdf



Important

Spectral calibration is required before running the test

To read the test results correctly, the sequencer needs to be calibrated with the DS-32 calibration kit (Applied Biosystems Cat. No. 4345831). This is a one-time calibration to set up correct spectral channels to read the test results. Refer to the DS-32 Matrix standards kit to prepare the DS-32 matrix standards. Run a Matrix Standard Set DS-32 (5FAM, JOE, NED, ROX) to perform a spectral calibration.

Thermal Cycling Programs:

Program 1 (PCR)

1 cycle	94°C 5 min
35 cycles	94°C 30 sec 52°C 30 sec 72°C 30 sec
1 cycle	72°C 5 min
	Hold at 4°C

Program 2 (PCR product clean up)

	37°C 25 min
	95°C 5 min
	Hold at 4°C

Program 3 (ST reaction)

1 cycle	94°C 4 minute
20 cycles	94°C 20 sec 55°C 30 sec 72°C 20 sec
	Hold at 4°C


Mutector™ II Assay Protocol:

A. PCR Amplification

Thaw all reagents and keep on ice. Spin down the reagents before use.

Mutation controls (CTL-BF) are required for each test run. A negative control (water) is suggested to run with the samples each time.

A.1. Collect 0.2 ml PCR strip tubes and label the tubes as follows:



Neg: Negative control (water)

Pos: Positive control (CTL-BF)

A.2. Prepare PCR reaction mix:

Master Mix = $18 \times \left(\frac{\text{Total number of sample}}{\text{Total number of sample}} + 2^* \right) \times 1.1^{**} = \text{_____ } \mu\text{L}$

PCR Primers = $1 \times \left(\frac{\text{Total number of sample}}{\text{Total number of sample}} + 2^* \right) \times 1.1^{**} = \text{_____ } \mu\text{L}$

*For negative and positive controls.

** For pipetting error.

Transfer entire volume of the reagents to one tube and mix the contents. This is the PCR reaction mix.

A.3. Add **19 μL** of PCR reaction mix into each tubes.

A.4. Add **1 μL** of nuclease-free water to the “**Neg**” tube.

A.5. Add **1 μL** of **CTL-BF** to the “Pos” tube.

A.6. Add **1 μL** of sample DNA (20-80 ng/μl) to each sample tube.

A.7. Place the PCR tubes in a thermal cycler and run **Program 1**.

<u>Program 1</u>	
1 cycle	94°C 5 min
35 cycles	94°C 30 sec
	52°C 30 sec
	72°C 30 sec
1 cycle	72°C 5 min
	Hold at 4°C



The procedure can be temporarily stopped after **Program 1**.
The PCR product can be stored at 4°C for next day test.

During the PCR amplification, prepare steps B1-B2.

B. PCR Product Clean-up

Prepare the C-UP mix

C-UP1 = 1.0 x number of samples x 1.1* = _____ μL

C-UP2 = 1.0 x number of samples x 1.1* = _____ μL

C-UP Buffer = 10 x number of samples x 1.1* = _____ μL

* For pipetting error.

B.1. Collect 0.2ml strip tubes, one tube for each PCR reaction.

Label the tubes the same as the PCR tubes.

B.2. Add 12 μL of **C-UP Mix** to each new tube.

B.3. Transfer 5 μL of PCR products to each corresponding tube.

B.4. Cap the tube, mix the contents and spin all tubes.

B.5. Incubate the tubes in a thermal cycler using **Program 2**.

<u>Program 2</u>
37°C for 25 min
95°C for 5 min
Hold at 4°C

During the incubation, prepare steps C1-C4.

C. ST Reaction (Mutation enrichment)

- C.1. Collect one 2ml tube and label the tube with “**ST**”. Prepare ST Mix using formula below:

ST Mix

$$\text{ST-BF} = 11 \times \left(\frac{\text{Total number of sample}}{\text{Total number of sample}} + 2^* \right) \times 1.1^{**} = \text{_____} \mu\text{L}$$


$$\text{DP-BF} = 2 \times \left(\frac{\text{Total number of sample}}{\text{Total number of sample}} + 2^* \right) \times 1.1^{**} = \text{_____} \mu\text{L}$$

*For negative and positive controls.

** For pipetting error.

Transfer entire volume of the reagents to one tube and mix the contents. This is the ST mix.

- C.2. Collect 0.2ml strip tubes, one for each PCR reaction. Label the tubes with sample ID. Label the tubes as follows:



- C.3. Add **13 μL** of **ST Mix** (from step C.1) into each of the tubes.
- C.4. Add **2 μL** of C-UP treated Negative PCR control to the “**Neg**” tube.
- C.5. Add **2 μL** of C-UP treated Pos to the “**Pos**” tube.
- C.6. Add **2 μL** of C-UP treated Sample to each corresponding sample tube.
- C.7. Mix the contents and spin all tubes.

- C.8.** Place the tubes into the thermal cycler and perform the ST reaction using **Program 3**.

Program 3

1 cycle	94°C 4 min
20 cycles	94°C 20 sec
	55°C 30 sec
	72°C 20 sec
	Hold at 4°C

D. ST Product Clean-up

Filter preparation

- D.1.** Collect the **TF - Filters** and **Collection Tubes** (one set for each ST reaction).
- D.2.** Snap off the bottom portion of the filter tip (ref. page 7 for snap off-line).
- D.3.** Centrifuge the TF - Filters at 1,000 x g (3000 rpm for most tabletop centrifuge) for 2-3 minutes to remove the excess liquid from the filters.
- D.4.** Discard the Collection Tubes and move the TF - Filters into a new Collection Tube. Label the Collection Tubes with sample ID. The TF - Filters are ready for use.
- D.5.** After the ST reaction, load all ST reaction contents (15µl) onto the top of the gel in each pre-prepared **TF - Filter**.
- D.6.** Centrifuge the **TF - Filters** at 1,000 x g (3000 rpm for most tabletop centrifuge) for 2-3 minutes.
- D.7.** Discard the **TF - Filters**. The solution in the tubes contains ST product and is ready for sample loading.

E. Sample Loading

E.1. Add **15 μ L** of the **Loading buffer** to each well of a sequencer adapter plate.

E.2. Transfer **2-4 μ L** of the filtered ST products into each well.

Signal may vary depending on the instrument used. It is recommended to adjust the loading volume to optimize the signal on your machine. If the signal is too strong, dilute the ST product with water (3-5 times) before loading the sample.

E.3. Load the plate to sequencer and run the pre-set Data Collection Program (ref. page 8).

F. Data Analysis

F1. Open the analysis software GeneMapper or GeneScan

F2. Follow the instructions to add the data for analysis.

The instructions are provided online:


GeneMapper: www.trimgen.com/docs/PartIII-Data-Analysis-GeneMapper.pdf

GeneScan: www.trimgen.com/docs/PartIV-Genescan.pdf.

F3. In the sample plot window (shows graphic data), zoom in the graphic by selecting area between size markers 25 (the 2nd marker) and 80 (the 6th marker) along the X-axis (see figure 1).

F4. Find the results for the mutation controls CTL-BF.

F5. The mutation controls CTL-BF will show a total of 5 peaks. Use these controls as a standard to identify peak(s) present in the samples.



F6. Sample analysis: A wild type (WT) peak is presented in every sample and serves as an internal control. Low peak height

indicates low concentration or poor quality of input DNA. If the peak is not observed, it indicates that the DNA amplification failed or the sample is 100% mutant, such as mutant cell lines.

When a sample contains a mutation, an additional peak shows up (sometimes there are multiple peaks due to an existence of multiple mutations). To identify the mutation, always compare the peak (size and color) with the peaks of the Mutation Controls (CTL-BF). Any peak that does not match with the mutant controls will not be considered.

The peak size may slightly shift due to migration differences between capillary tubes. The size shift can be identified by comparing the wild type peak in the sample to the wild type peak in the Mutation Controls CTL-BF. When the shift occurs, the peak size for the mutations will have a proportional shift.

G. Troubleshooting

G1. “Color leak-through”

When the sample DNA concentration is too high, the ST reaction generates a strong fluorescent signal >5,000 rfu. Fluorescence spillover will occur. For example, the black peak of the wild type signal may be observed in the red and/or blue channels. This color spillover is caused by limitation of the instrument. The “leak-through” peak will have the exact same peak size as the original peak. Because the mutation peaks have different peak size, leak-through will not affect data analysis.

G2. The peak signal is too high

The assay is set at a condition to detect mutations in a small sample, such as DNA extracted from fine needle aspiration (FNA) sample. For regular FFPE sample, the assay signal may be too high to analyze (peak height >8000 rfu, can not see the top of the peak or the peak is highlighted with pink color). Diluting the final STA product with de-ionized water can efficiently reduce the signal and optimize the peak height. Do not dilute the assay reagents, it will cause improper enzymatic reaction and generate a miss call. Each laboratory has different PCR instrument(s), the signal intensity may vary among the laboratories, first time users should perform a titration of final STA product. Once the dilution factor is determined, the assay will show consistent results.

G3. Graphic data will not automatically show

Check the raw data. If the signals from the sample and size standards are too low, the capillary tube may be blocked by a bubble. The sample needs to be re-loaded. The ST products will compete with the size standard DNA to enter the capillary tube.

The high sample signal will reduce the size standard signal. Diluting the final ST product with de-ionized water will easily resolve this problem.

The size standard may be miscalculated. Check the size standard and manually correct the size standard and reanalyze the data.

G4. Peak size shift

Sometimes the peak size may vary between the capillary tubes. Air bubbles or capillary tube conditions could cause the shift. Comparing the wild type peaks in the sample and Mutation Controls can identify the peak size shift. If the shift occurs, the mutation peaks will show a similar shift.

G5. Extra peaks do not match any peaks in CTL-BF controls.

In some circumstances, for example sample DNA is too low or PCR is not amplified properly, the PCR primers may form primer dimers. These primer dimers can generate blue peak(s) with different sizes. Any peaks that do not match the CTL-BF controls are not considered a mutation signal.

G6. Background noise

The background of the test is normally low and high background noise may be caused by poor quality DNA. When the PCR reaction does not perform properly, the intensity of the wild type peak is usually lower than 500 and many small peaks are observed. These peaks are not considered as signal.

G7. Mutation peak cut-off

For some samples, a small peak may be observed in one of the mutation positions. To verify the peak, you need to confirm the signal strength of the wild type peak. If the wild type peak is too high (can not see the top of the peak and the peak is highlighted with pink color), your ST reaction is too strong and the small peak may be "pull up" from background noise. Follow trouble shooting 2. to dilute the final product of the ST reaction with de-ionized water. After dilution, reload the sample. If you can see the top of the wild type peak, use the following calculation to identify the small peak:

$$\text{Ratio} = (\text{Area of mutant peak}) / (\text{Area of wild type peak})$$

If the ratio is larger than 0.06, the peak is determined to be a mutation peak (the ratio does not represent the percentage of the mutation in the sample). Otherwise, the peak is a background pull-up and does not indicate the presence of a mutation in the sample.