

# Swab and Go...

*a quick method for DNA sample collection*

## SwabPack™

**Buccal cell collection and shipping kit**

The kit includes non-invasive sample collection swabs, and sample shipping envelope. Easy-Swab™ provided in the kit is specifically designed to efficiently collect buccal cells. The DNA recovery from the Easy-Swab™ is much higher than with traditional swabs. The swab can be used to collect many types of cells such as buccal or nasal cells. After collection, the swab can be dried at room temperature and then shipped to a laboratory for testing. The swab holder and shipping envelope are designed for shipping the swab samples using standard mail so no special shipping procedure is required.

## SwabPack™



**1-888-825-6005**  
**www.trimgen.com**

34 Loveton Circle, Suite 210 Sparks, MD 21152  
1-8-34 Oyodonaka Lota-Ku Osaka, Japan

**Easy to Collect**  
**Easy to Ship**  
**Easy to Recover DNA**

### Maximal cell collection

The Easy-Swab™ is designed with special materials and a unique surface that can collect cells with high efficiency.

### Long term storage is possible

After collection, the Easy-Swab™ can be dried and keep at room temperature for up to several months before the DNA or RNA extraction is performed.

### Save sample shipping cost

The collected sample can be shipped at room temperature using standard mail. No special shipping procedure is required.

### Recover DNA in 5 minutes

The DNA can be recovered from the Easy-Swab™ in 5 minutes using TrimGen's BuccalQuick™ DNA kit (Cat No. BQ-100). The DNA extracted from a single swab is sufficient for performing more than 50 PCR reactions.

### Product Information

Catalog Number	Kit Size
SP-50	50 samples
SP-100	100 samples

For Research Use Only.  
Not for use in diagnostic procedures.