

Mutector™ Reagents



CYP2C19

Genotyping



Multiplex

Multiplexing

Sequencing accuracy

Results in less than 4 hours

Target SNPs

SNP	Variation	SNP	Variation
CYP2C19 * 2	19154G > A	CYP2C19 * 6	12748G > A
CYP2C19 * 3	17948G > A	CYP2C19 * 7	19294T > A
CYP2C19 * 4	1A > G	CYP2C19 * 8	12711T > C
CYP2C19 * 5	90033C > T	CYP2D6 * 17	-806C > T



FEATURES

Clear-Cut Results

- Individual SNP is defined by peak's color & size
- Easy to analysis the allele

Multiplexing Capability

- All SNPs in different exons can be tested in two tubes

Flexibility

- Easy to add new SNPs without changing protocol

Simple Protocol for Same Day Results

- Have results in 3-4 hours

50% off

**Validation
Kit**

ORDER INFORMATION

Product Name	Cat. No.	Kit Description	Size
CYP2C19 Genotyping Kit	GP12	Mutector™ Reagents for CYP2C19 Genotyping	32 rxn

For research use only, not for use in diagnostic procedures

888.825.6005

www.trimgen.com

34 LOVETON CIRCLE, SUITE 210, SPARKS, MD 21152 USA,
PHONE: 410.472.1100 FAX: 410-472-1303